



Amidite Diluent Station

We offer a station to dilute amidites. The station is used to dilute a single amidite at a time. We can supply the station with any bottle cap for the amidite. Some examples are:

- Clip Cap used on ABI instruments to hold septum top bottles
- Screw cap, 20 mm used on Expedites
- Screw Cap, 24 mm used on some High Throughput Synthesizers
- Screw Cap, 28 mm used on some High Throughput Synthesizers

If you need a cap that is not listed above, let us know and chances are that we can supply it.


The cap for the Acetonitrile Bottle can be any of the following. You will need to specify which one you want:

- A 38 mm cap used on 2 & 4 liter bottles
- A 45 mm cap used on the 2.5 liter bottles used in Europe.
- A 28 mm cap for a 450 ml bottle

We can make this product in a lot of different variations. If you need something other than our standard off the shelf version. Some possible variations are

- More than one amidite bottle fill at the same time
- Different Amidite Bottle caps on the same Diluent Station (fill one at a time but can have 2 or 3 different type of amidite bottle caps).

Part Number is for our Off the Shelf Amidite Diluent Station: 404-101687





Shown with 28 mm
screw cap for amidite
bottle

Optional Clip Cap for Amidite
Bottle.